

Ultrplex® 360


B | BRAUN

LLDorder 7922 - lw 10289197 (vormals gw 230424)

Peripheral
Spinal
Epidural

en	Ultrasound-visible needle for peripheral nerve blocks	4
de	Ultraschallsichtbare Kanüle für periphere Nervenblockaden	5
bg	Видима под ултразвук канюла за периферни нервни блокади	7
cs	Jehla pro blokády periferních nervů viditelná pod ultrazvukem	8
da	Ultralydsynlig nål til perifere nerveblokader	10
el	Βελόνα για αποκλεισμό των περιφερικών νεύρων, κατάλληλη για χρήση με υπέρηχο	11
es	Cánula visible mediante ecografía para bloqueos nerviosos periféricos	13
et	Ultraheliks sobiv nõel perifeersetele närviblokaadidele	14
fi	Ultraäänikäyttöön sopiva neula ääreishermosten johtopuudutukseen	16
fr	Aiguille visible sous échographie pour blocs nerveux périphériques	17
hu	Ultrahangot visszaverő tű perifériás idegblokkoláshoz	19
it	Ago ecogenico per blocchi dei nervi periferici	20
lt	Ultragarsu kontroliuojama adata periferinei nervų blokada	22
lv	Ar ultraskaņu redzama adata perifērā nerva blokādei	23
nl	Bij echografie zichtbare naald voor perifere zenuwblokkades	25

no	Ultralydegned nål for perifere nerveblokkeringer	26
pl	Igła do blokad nerwów obwodowych o wyższej echogeniczności przeznaczona do stosowania z USG	28
pt	Agulha visível por ultra-sons para bloqueio de nervos periféricos	29
pt	BR Agulha visível por ultrassom para bloqueio de nervos periféricos	31
ro	Ac cu vizibilitate ecografică pentru stimularea blocuri nervoase periferice	32
sk	Ultrazvuková ihla na blokády periférnych nervov	34
sl	Igla za izvajanje periferne živčne blokade, primerna za ultrazvočno opazovanje	35
sv	Med ultraljudsynlig nål för perifera nervblockeringar	37
th	เข็มสำหรับระงับการทำงานของประสาทรอบนอกชนิดมองเห็นด้วยอัลตราซาวด์ได้	38
tr	Periferik sinir bloğu için ultrasonda görülebilen iğne	40
zh	用于外周神经阻滞的超声可视针头	41

en Instructions for Use

Device description

Luer-connector

The female luer connector of the Ultraplex® 360 devices is a standardized connector. It can be connected to syringes with male luer connector. The syringes can be used for flushing the needle, aspirating and infusing anesthetics.

Assessment of Magnetic Resonance (MR) compatibility

Ultraplex® 360 Needle: The Device has not been evaluated for safety in the MR environment, since it is not expected to enter the MR Environment in a typical application. Performing an MR exam on a person who has this medical device may result in injury or device malfunction.

For further technical characteristics such as needle size, see product packaging.

Sterility

Sterilized by EO.

Do not use if package is damaged.

Visually inspect the sterile barrier system for breaches and integrity prior to use.

Intended purpose

Ultrasound-visible needle for peripheral nerve blocks.

Indication

Plexus and peripheral nerve blocks using ultrasound respectively for precise detection of definite nerve areas. Nerve blocks are used as anesthesia for surgical procedures as well as for pain therapy involving the upper and lower extremities with various block techniques.

Patient population

All patients in need for anesthesia or pain therapy by plexus or peripheral nerve blocks. If, despite adequate explanation, the patient refuses plexus anesthesia, their wishes should be respected. Likewise, mentally handicapped patients and those with psychiatric problems need careful preoperative assessment.

Intended User

The product should only be used by healthcare workers that have been adequately trained in this technique.

Contraindications

- Known pre-existing neurological diseases
- Infection / inflammation at the desired site of insertion / puncturing
- Excessive preoperative anxiety (not responsive to interventions) and inability to tolerate positioning
- Pre-existing coagulopathies, either endogenous (e.g. thrombocytopenia) or iatrogenic (e.g. warfarin therapy) due to increased risk of hematoma
- Anatomic anomalies that make needle placement difficult via ultrasound and anatomical landmarks

- Hepatic disease that may interfere with clearance
- Failure to obtain patient's consent
- Known hypersensitivity to one or more materials used
- Contraindications for axillary blocks: lymphangitis and lymphadenectomies
- For information on other contraindications against providing regional anesthesia, please refer to the specialist literature.

Residual Risks/Side effects

- Nervous lesions
- Pneumothorax
- Pain and hematoma at the site of insertion of the needle
- Incomplete block
- Local anesthetics related toxicity: cardiac, neurologic, allergic
- Incorrect epidural or intrathecal injection of local anesthetic in the case of peripheral blocks in the area of the spinal column
- Infection (sepsis)

Warning

- Single use only. Do not re-use.
- Do not resterilize.
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- If paresthesia is caused inadvertently as a result of direct contact with the nerve, the needle must not, under no circumstance, be inserted any further.
- An unusually high injection pressure may indicate that an injection has been administered intraneurally. Discontinue the injection, otherwise irreversible nerve damage may occur.
- Visually inspect the sterile barrier system for breaches and integrity prior to use.
- Do not use if package is damaged.
- Use the properly stored product by the date stated on the packaging.

Operating Instructions

Use strictly under aseptic conditions.

Disinfect block area thoroughly. Repeat this at least two times. Allow time to act between each application and cover with a sterile fenestrated drape.

Position the patient and inject a skin wheal of local anesthetic at the puncture site and if necessary, infiltrate the puncture track.

Prime tubing and needle with solution.

Advance the needle under ultrasound vision towards the nerve while continually collaborating with the patient and their feedback on any paralgesia or paresthesia. Test injections with 0.9% isotonic NaCl support verification of the needle position. Continue this procedure until the nerve block needle is situated in appropriate distance to the nerve. Distinct reflexes of the tip

geometry and the specific pattern of the distal shaft area of the needle should be visible and indicate the needle position in relation to the nerve. Please note that the visibility may depend on a number of factors such as proper setting of the ultrasound unit, patient related conditions of the tissue, artefacts and proper alignment of the US probe and the needle.

Aspirate for possible intravascular placement.

Following negative aspiration, a test dose of local anesthetic can be injected through the extension tubing.

Duration of use

Sterile device for single use.

The Ultraplex® nerve block needle must only be used for detecting the target nerve and for injecting the local anesthetic. It must then be removed.

Disposal

Disposal according to local guidelines and clinical protocols.

Storage and handling conditions

Please refer to information on packaging.

Keep away from sunlight.

Keep dry.

Notice to the user

If, during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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AU For sterilization type see primary packaging
B. Braun Australia Pty Ltd
Level 2, 4 Burbank Place
Norwest NSW 2153
Australia

CA Distributed by: **B. Braun Medical Inc.**
Bethlehem, PA USA 18018-3524

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US Rx only
Pictorials are for reference only
Distributed by: **B. Braun Medical Inc.**
Bethlehem, PA 18018-3524
Not made with DEHP

Manufacturing site:
B. Braun Aesculap Japan Co., Ltd.
Tochigi Factory - Hospital Care
285 Ogaki, Tsuga-machi, Tochigi-shi,
Tochigi, 328-0101 Japan

0123

0126 X0620519

B. Braun Melsungen AG
34209 Melsungen
Germany
www.bbraun.com

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B.Braun (Thailand) Ltd.,