

Absorbent underpad with SAP, non-sterile

TD-03-I.2.c-1.5

Revision 1

date: 2024-07-08

amendment of: N/A

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Po

ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com

CLASSIFICATION

Class I medical device compliant with Medical Devices Regulation (EU) 2017/745 (MDR)

INTENDED PURPOSE

 disposable device is used to absorb body fluids that contain potential pathogens that pose a risk of infection to the patient. Ethylene oxide sterilization is required for sterile use.

FEATURES

- quilted absorbent insert
- composed of layers: polypropylene non-woven, cellulose wadding, cellulose pulp with superabsorbent, cellulose wadding, PE foil
- white absorbent insert
- blue foil layer
- single-use
- non-sterile
- shelf life: 5 years
- packaging sealed plastic bag



CATALOGUE NUMBER	DIMENSIONS	COLOUR	TYPE OF PACKAGING	SALES PACKAGING	Transport Packaging
BT-PCS4060	40 x 60 cm	white/blue	foil	25 pcs	4 x 25 pcs
BT-PCS6060	60 x 60 cm	white/blue	foil	25 pcs	4 x 25 pcs
BT-PCS6090	60 x 90 cm	white/blue	foil	25 pcs	4 x 25 pcs