

CARESITE Luer Access Device

Making access clear and easy for improved patient care.



IV Accessories

CARESITE Luer Access Device

Helppo tarkistaa ulkoisesti



Helppo puhdistaa.



Hyvä tartuntapinta.



Helppo yhdistää.



CARESITE tuoteominaisuudet

Virtausnopeus	208 ml/min
Paine	300 psi @ 10mL/sec
Jäännösvolyymi	0.20 mL
Täyttövolyymi	0.22 mL
Nesteen siirto	0.03 mL (positive)
Lävistettävyyys	216 krt/min
✓ Suljettu systeemi	✓ Lateksiton
✓ Neulaton	✓ Veri / lipidi yhteensopiva
✓ Split septum teknologia	✓ Luer lock / Luer slip yhteensopiva

Edut käyttäjälle:

Suojaa mikrobiologiselta kontaminaatiolta

- Sileä, pyyhittävä / puhdistettava
- Selkeä muotoilu
- Läpinäkyvä

Estää katetrin tukkeutumisen

- Positiivinen venttiili (venttiilin rakenne)

Helppokäyttöinen

- Korkeapaineenkestävä
- Muotoilu edistää helppokäyttöisyyttä
- Lääkeyhteensopiva

Description	Length (cm)	Priming volume (ml)	Tubing Inner ø (mm)	Latex free	DEHP free	Sales Unit / Pcs.	Code No. (REF)
Caresite 	-	0.22	-	▪	▪	200 / 100**	415122-01
Extension Sets							
Small bore Extension set with Caresite, Spinlock (high pressure – for power injectors) 	20	0.5	1.3	▪	▪	100	470100-01
Standard bore Extension set with Caresite, Spinlock 	15	0.9	2.8	▪	▪	100	
Bifurcated small bore Extension set with 2 Caresite, Spinlock 	18	0.9	1.3	▪	▪	50	470106-01

* lipids and cremophor
** dispenser box

Detailed information for healthcare professionals

For more information about risk prevention in infusion therapy, visit:

www.safeinfusiontherapy.com