


# Sol-M™ Blunt Fill Needle


Sol-M™ Blunt Fill Needles were designed to reduce the risk of needle-stick injury by removing sharps while drawing up and preparing medication. For use when drawing medication and can also be used for adding medication in needleless IV systems.



CE 2797

## Features:

- Eliminates the unnecessary use of sharps
- Color-coded hubs for easy identification
- Needle hub is compatible with all Luer Lock and Luer Slip Tip syringes
- Available with 5-micron filter designed to restrain drawing of shattered particles when using ampules (REF: 110022F)
- Not manufactured with Latex, PVC and DEHP



# Sol-M™ Blunt Fill Needle



5-micron filter needle may be used with glass ampoules.

Designed to prevent even 5-micron particles from entering the syringe.

Note: In case particle are visible in the vial or ampoule it is recommended to discard the medication to prevent injections of those particles in the patient.

## Sol-M™ Blunt Fill Needles

Ref. Number	Product Description	Gauge	Length (inches)	Length (mm)	Box Qty	Case Qty
110021	Sol-M Blunt Fill Needle 18G x 1"	18G	1"	25 mm	100	1000
110022	Sol-M Blunt Fill Needle 18G x 1 1/2"	18G	1 1/2"	40 mm	100	1000
110023	Sol-M Blunt Fill Needle 18G x 2"	18G	2"	50 mm	100	1000
110022F	Sol-M Blunt Fill Needle w/Filter 18G x 1 1/2"	18G	1 1/2"	40 mm	100	1000
110023F	Sol-M Blunt Fill Needle w/Filter 18G x 2"	18 G	2"	50 mm	100	1000

Please contact your local Sol-Millennium® representative for additional sizes of the Sol-M™ Blunt Fill Needle.

**Intended use:** Sol-M™ Blunt Fill Needle is used in conjunction with a syringe as an additive device for aspiration from multi-dose medicine vials or injection into I.V. Systems and pre-slit septum covering injection sites.

**Intended user:** Licensed healthcare professionals (HCP).

**Cautions:**

- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/ injury.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.

This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN.

Sol-Millennium Europe Sp. z o.o. All rights reserved.

EMA-MD-437-0323

**For more information, please contact:**



Sol-Millennium Europe Sp. z o.o.  
Warsaw HUB, Rondo Daszyńskiego 2B  
00-843 Warsaw, Poland



info.eu@sol-m.com



**Sol M**

Building  
a Healthier  
Tomorrow