

Sol-Care™ Safety Blood Collection Needle

- Sol-Care™ Safety Blood Collection Needle is a safety engineered device. The safety mechanism can be activated immediately after blood draw, helping to protect you against needlestick injury.



Sliding safety shield

Available with
or without
pre-attached
holder

Blood flashback
and flow are
easy to visualize

CE 0123

Sol-Care™ Safety Blood Collection Needle

Sol-Care™ Safety Blood Collection Needle with Luer Adapter

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (inches)	Needle Length (mm)	Tubing Length (inches)	Tubing Length (mm)	Box Qty	Case Qty
110201030035	25 G	0.5 mm	3/4"	19 mm	12"	305 mm	50	200
110201030033	23 G	0.6 mm	3/4"	19 mm	12"	305 mm	50	200
110201030031	21 G	0.8 mm	3/4"	19 mm	12"	305 mm	50	200
110201030034	25 G	0.5 mm	3/4"	19 mm	7"	178 mm	50	200
110201030032	23 G	0.6 mm	3/4"	19 mm	7"	178 mm	50	200
110201030030	21 G	0.8 mm	3/4"	19 mm	7"	178 mm	50	200

Sol-Care™ Safety Blood Collection Needle with pre-attached holder

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (inches)	Needle Length (mm)	Tubing Length (inches)	Tubing Length (mm)	Box Qty	Case Qty
110201030010	25 G	0.5 mm	3/4"	19 mm	12"	305 mm	25	200
110201030006	23 G	0.6 mm	3/4"	19 mm	12"	305 mm	25	200
110201030002	21 G	0.8 mm	3/4"	19 mm	12"	305 mm	25	200
110201030012	25 G	0.5 mm	3/4"	19 mm	7"	178 mm	25	200
110201030008	23 G	0.6 mm	3/4"	19 mm	7"	178 mm	25	200
110201030004	21 G	0.8 mm	3/4"	19 mm	7"	178 mm	25	200

Intended use: Used for blood collection or/and short-term infusion of intravenous fluids (when used without the Holder with Male Luer Adapter accessory).

Used for blood collection (when used with the Holder with Male Luer Adapter accessory).

Intended user: Licensed healthcare professionals (HCP).

Cautions:

- Do not use the device within the Magnetic Resonance (MR) Environment.
- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/ injury.
- Sharps object. Risk of sharps injuries.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.

This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN.

Sol-Millennium Europe Sp. z o.o. All rights reserved.

EMA-BC-611-1123

For more information, please contact:



Sol-Millennium Europe Sp. z o.o.
Warsaw HUB, Rondo Daszyńskiego 2B
00-843 Warsaw, Poland



info.eu@sol-m.com



Sol M Building a Healthier Tomorrow