

elastoDERM

foil dressing, self-adhesive, sterile

features:

- made of polyurethane film
- adhesive
- covered with hypoallergenic acrylic glue
- transparent - allows observation of the dressed place without removing the dressing
- vapour permeability guarantees correct gas exchange between the foil and skin and prevents it from excessive accumulation of moisture and maceration
- barrier against microorganisms and bacteria minimizing infection risk
- very thin and elastic - adjusts to the body shape and guarantees freedom of movement
- may remain on the skin up to 7 days; optimum holding time of the dressing is from 3 to 4 days
- the application system - "frame" or "three-stage" type - facilitates precise and painless application of the dressing
- available in various sizes and options:
 - with central absorbent pad
 - with incision which allows adjusting of the dressing to cannula
 - without incision or absorbent pad
- waterproof
- hypoallergenic
- sterile



elastoDERM F

dressing without absorbent pad with "frame" application type

ref.	dimensions	shape	indication	packaging unit/ pouch	intermediate packaging/box	bulk packaging
812003	4,4 cm x 4,4 cm		- secondary dressing - bedsore prevention, works as the so called "second skin" - for small epidermal abrasions and damages on the skin surface - fixing and protecting cannulas inserted into central veins	1 pc.	100 pcs	36 boxes
812005	6 cm x 7 cm		- secondary dressing - bedsore prevention, works as the so called "second skin" - for small epidermal abrasions and damages on the skin surface - fixing and protecting cannulas inserted into central veins	1 pc.	100 pcs	30 boxes
with description field						
812008	10 cm x 12 cm		- secondary dressing - bedsore prevention, works as the so called "second skin" - for small epidermal abrasions and damages on the skin surface - fixing and protecting central veins catheters - fixing epidural catheters	1 pc.	100 pcs	20 boxes
812011	10 cm x 15 cm		- secondary dressing - bedsore prevention, works as the so called "second skin" - for small epidermal abrasions and damages on the skin surface - fixing and protecting central veins catheters - fixing epidural catheters	1 pc.	100 pcs	20 boxes
with description field						

elastoDERM F-IV

IV cannula dressing with „frame" application type

812006	6 cm x 7 cm		- fixing and protecting peripheral venous cannulation in adults - protecting cannulas inserted into central veins	1 pc.	100 pcs	30 boxes
with description field						

elastoDERM

dressing without absorbent pad with "three-stage" application type

811210	6 cm x 7 cm		- secondary dressing - bedsore prevention, works as the so called "second skin" - for small epidermal abrasions and damages on the skin surface - fixing and protecting cannulas inserted into central veins	1 pc.	100 pcs	30 boxes
--------	-------------	--	---	-------	---------	----------

elastoDERM PAD

dressing with central absorbent pad with "three-stage" application type

802012	6 cm x 7 cm		- protection of postoperative wounds and stitches	1 pc.	100 psc	12 boxes
802014	9 cm x 10 cm			1 pc.	50 pcs	24 boxes
802015	9 cm x 15 cm			1 pc.	50 pcs	24 boxes
802016	10 cm x 20 cm			1 pc.	50 pcs	18 boxes
802018	10 cm x 25 cm			1 pc.	50 pcs	15 boxes