

# Sol-Care™ Safety Blood Collection Needle

- Sol-Care™ Safety Blood Collection Needle is a safety engineered device. The safety mechanism can be activated immediately after blood draw, helping to protect you against needlestick injury.



Sliding safety shield

Available with  
or without  
pre-attached  
holder

Blood flashback  
and flow are  
easy to visualize

CE 0123

# Sol-Care™ Safety Blood Collection Needle

## Sol-Care™ Safety Blood Collection Needle with Luer Adapter

| Ref. Number  | Gauge | Needle outer diameter (mm) | Needle Length (inches) | Needle Length (mm) | Tubing Length (inches) | Tubing Length (mm) | Box Qty | Case Qty |
|--------------|-------|----------------------------|------------------------|--------------------|------------------------|--------------------|---------|----------|
| 110201030035 | 25 G  | 0.5 mm                     | 3/4"                   | 19 mm              | 12"                    | 305 mm             | 50      | 200      |
| 110201030033 | 23 G  | 0.6 mm                     | 3/4"                   | 19 mm              | 12"                    | 305 mm             | 50      | 200      |
| 110201030031 | 21 G  | 0.8 mm                     | 3/4"                   | 19 mm              | 12"                    | 305 mm             | 50      | 200      |
| 110201030034 | 25 G  | 0.5 mm                     | 3/4"                   | 19 mm              | 7"                     | 178 mm             | 50      | 200      |
| 110201030032 | 23 G  | 0.6 mm                     | 3/4"                   | 19 mm              | 7"                     | 178 mm             | 50      | 200      |
| 110201030030 | 21 G  | 0.8 mm                     | 3/4"                   | 19 mm              | 7"                     | 178 mm             | 50      | 200      |

## Sol-Care™ Safety Blood Collection Needle with pre-attached holder

| Ref. Number  | Gauge | Needle outer diameter (mm) | Needle Length (inches) | Needle Length (mm) | Tubing Length (inches) | Tubing Length (mm) | Box Qty | Case Qty |
|--------------|-------|----------------------------|------------------------|--------------------|------------------------|--------------------|---------|----------|
| 110201030010 | 25 G  | 0.5 mm                     | 3/4"                   | 19 mm              | 12"                    | 305 mm             | 25      | 200      |
| 110201030006 | 23 G  | 0.6 mm                     | 3/4"                   | 19 mm              | 12"                    | 305 mm             | 25      | 200      |
| 110201030002 | 21 G  | 0.8 mm                     | 3/4"                   | 19 mm              | 12"                    | 305 mm             | 25      | 200      |
| 110201030012 | 25 G  | 0.5 mm                     | 3/4"                   | 19 mm              | 7"                     | 178 mm             | 25      | 200      |
| 110201030008 | 23 G  | 0.6 mm                     | 3/4"                   | 19 mm              | 7"                     | 178 mm             | 25      | 200      |
| 110201030004 | 21 G  | 0.8 mm                     | 3/4"                   | 19 mm              | 7"                     | 178 mm             | 25      | 200      |

**Intended use:** Used for blood collection or/and short-term infusion of intravenous fluids (when used without the Holder with Male Luer Adapter accessory).

Used for blood collection (when used with the Holder with Male Luer Adapter accessory).

**Intended user:** Licensed healthcare professionals (HCP).

### Cautions:

- Do not use the device within the Magnetic Resonance (MR) Environment.
- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/ injury.
- Sharps object. Risk of sharps injuries.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.

This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN.

Sol-Millennium Europe Sp. z o.o. All rights reserved.

EMA-BC-611-1123

### For more information, please contact:



Sol-Millennium Europe Sp. z o.o.  
Warsaw HUB, Rondo Daszyńskiego 2B  
00-843 Warsaw, Poland



info.eu@sol-m.com



**Sol M** Building a Healthier Tomorrow