

# Sol-Care™ Safety Multi-sample Needle

- Safety-engineered tool for multi-sample venous blood collection. The device is designed with a safety mechanism covering the needle after use. In the activated position, the needle cover guards against accidental needlestick during normal handling and disposal. It is also offered with a pre-attached holder.



CE 0123

# Sol-Care™ Safety Multi-sample Needle

## Sol-Care™ Safety Multi-sample Needle

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (inches)	Needle Length (mm)	Box Qty	Case Qty
SMSN22125	22 G	0.7 mm	1 1/4"	32 mm	48	480
SMSN2215	22 G	0.7 mm	1 1/2"	38 mm	48	480
SMSN21125	21 G	0.8 mm	1 1/4"	32 mm	48	480
SMSN2115	21 G	0.8 mm	1 1/2"	38 mm	48	480
SMSN20125	20 G	0.9 mm	1 1/4"	32 mm	48	480
SMSN2015	20 G	0.9 mm	1 1/2"	38 mm	48	480



## Sol-Care™ Safety Multi-sample Needle with Pre-attached Holder

Ref. Number	Gauge	Needle outer diameter (mm)	Needle Length (inches)	Needle Length (mm)	Box Qty	Case Qty
SMSNH22125	22 G	0.7 mm	1 1/4"	32 mm	50	400
SMSNH2215	22 G	0.7 mm	1 1/2"	38 mm	50	400
SMSNH21125	21 G	0.8 mm	1 1/4"	32 mm	50	400
SMSNH2115	21 G	0.8 mm	1 1/2"	38 mm	50	400
SMSNH20125	20 G	0.9 mm	1 1/4"	32 mm	50	400
SMSNH2015	20 G	0.9 mm	1 1/2"	38 mm	50	400



**Intended use:** Used with tube holder as a system in routine venipuncture/ blood collection procedure. In the activated position, the safety shield protects against accidental needle stick injury during normal handling and disposal.

**Intended user:** Licensed healthcare professionals (HCP).

### Cautions:

- Do not forcefully de-activate or re-activate the safety mechanism after it has been activated.
- Single use device. Re-use or use if the package is damaged may lead to infection or other illness/injury.
- Sharps object. Risk of sharps injuries.

The content of this material is directed at and restricted to a professional audience only (Healthcare Professionals). This material is not medical advice, guideline or a recommendation. Before using any particular product always read the package insert, product label and/or instructions for use.

This material is distributed by Sol-Millennium Europe Sp z o.o. with a registered address at Twarda 18, 00-105 Warsaw, Poland, Registry Court: District Court for Warsaw, XIIth Economic Department KRS 0000434783, EU VAT: PL5252535720, Share capital: 500.000 PLN.

Sol-Millennium Europe Sp. z o.o. All rights reserved.

EMA-BC-612-1123

### For more information, please contact:



Sol-Millennium Europe Sp. z o.o.  
Warsaw HUB, Rondo Daszyńskiego 2B  
00-843 Warsaw, Poland



info.eu@sol-m.com



**Sol M** Building  
a Healthier  
Tomorrow