

# Veterinary adhesive tape with *fun* paw prints!



## PawTape

### VETERINARY ADHESIVE TAPE

- Strong - 33 lbs. tensile strength
- Synthetic adhesive tape
- Crisp, easy hand-tear
- Delivers consistent support
- Easy to apply
- Made in the USA

NEW!



### Product Applications

- Immobilize wounded areas
- Splinting, strapping, and patient positioning
- Hold other dressings in place
- Protection and support

### PawTape

Ordering Information...

Catalog #	Size	Rolls/Box	Boxes/Case	Rolls/Case
291-005-288	½" x 10 yds.	24	12	288
291-010-144	1" x 10 yds.	12	12	144
291-020-072	2" x 10 yds.	6	12	72

### How to Use:



After cleaning and bandaging the area of application, place the adhesive tape half on the bandage and half on the skin.



Apply as necessary.



Hand-Tear or cut with scissors and place hand over tape for upwards of 30 seconds to slightly warm the adhesive.