



2-LAYER SURGICAL DRAPE
with central fenestration
sterile
TD-11

att.no.3 pt.2 subpt.2.3

revision 2

date: 2017-10-19

amendment of: 2023-02-01

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.

ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- Patient and equipment drapes for surgical procedures
- Intended for use in the operating room
- Intended for use during surgical procedures
- As an impermeable barrier for liquids and microorganisms
- Absorb excess fluids during surgical procedures
- Elimination of dusting due to the use of cotton-free materials

FEATURES:

- Two-layer surgical drape ALPHAtex, made of 40 g/m² laminate
- The non-woven layer absorbs moisture and fluids, the PE film layer prevents strike-through of fluids
- Sterile, single-use
- Providing an optimal barrier against microbial migration
- Fulfill the requirements of the standard EN 13795-1
- The fenestration facilitates the application of the drape, allows a more accurate fit of the drape to the surgical site, which increases the effectiveness of the operating field protection and ergonomics of the work of medical personnel
- Non-adhesive fenestration
- Paper-foil packaging
- Label with four adhesive tag labels for data archiving. The following information is available on each adhesive tag label: ref number, expiry date, batch number, manufacturer's data and barcode



SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	FENESTRATION	INDIVIDUAL PACKAGING	SALES PACKAGING	TRANSPORT PACKAGING
AT-NFF-S 3 E	45cm x 60cm	∅ 6cm	1 pc.	100 x 1 pc.	4 x 100 pcs
AT-NFF-S 2 E	45cm x 75cm	6cm x 8cm	1 pc.	100 x 1 pc.	4 x 100 pcs
AT-NFF-S 4 E	75cm x 90cm	∅ 7cm	1 pc.	50 x 1 pc.	4 x 50 pcs
AT-NFF-S 5 E	75cm x 90cm	6cm x 8cm	1 pc.	50 x 1 pc.	4 x 50 pcs