

PRODUCT INFORMATION


VACUETTE® Blood Collection Tubes

ESR Sodium Citrate

FILL VOLUME CHART

Tube type:

VACUETTE® ESR TUBES GLASS



The fill volume chart assists in a visual control of the specimen volume. This chart is prepared on a 1:1 scale. Please note that some printers are not able to print the correct scale.

The fill tolerance of the **VACUETTE®** ESR Tubes is within $\pm 10\%$ of the stated volume according to the international standards ISO 6710 and CLSI.