



paraffiNET
paraffin gauze dressing, sterile

TD-14

att.no.3 pt.2 subpt.2.3

revision 2

date: 2016-06-03

amendment of: 2023-02-01

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.


ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02,06 fax. +48 32 370 38 08

INDICATION:

- keeping proper moisture in the wound
- preventing the wound edges from drying and maceration
- recommended for uninfected wounds
- dressing chronic wounds: bedsores, ulcers, etc.
- treating surgical wounds, burn wounds, places of injections, skin transplants and damaged skin
- protection against postradiation skin damage

FEATURES:

- wide-mesh dressing made of cotton gauze covered with paraffin
- structure of the carrier material allows exudate to be freely drained to the covering secondary dressing
- covering with a layer of paraffin reduces risk of adhesion to the wound and allows a painless change of dressing
- neutral properties of paraffin allow for using the dressing in combination with other preparations containing medicinal substances or antiseptics
- can be shaped or cut according to the needs aseptically
- needs to be covered with a dressing with absorbent properties, e.g. compress, high absorption dressing, polyurethane foam
- non-adhesive, necessary fixing with a tape, bandage or a net tubular bandage
- hypoallergenic
- sterile
- shelf life: 5 years



SIZE AND REFERENCE

CATALOGUE NUMBER	DIMENSIONS	INDIVIDUAL PACKAGING PAPER-ALUMINUM FOIL	SALES PACKAGING BOX	TRANSPORT PACKAGING
PAR-01	5 cm x 5 cm	1 pc.	10 pcs	50 x 10 pcs
PAR-02	10 cm x 10 cm	1 pc.	10 pcs	50 x 10 pcs
PAR-06	10 cm x 40 cm	1 pc.	10 pcs	50 x 10 pcs
PAR-07	15 cm x 20 cm	1 pc.	10 pcs	50 x 10 pcs