

NONVI lux

non-woven swab, non-sterile

features:

- made of hydrophilic medical non-woven material, 30 or 40 g/m²
- 4 plies
- soft and non-fraying
- can be sterilized with steam, ethylene oxide or gamma radiation

indication:

- after sterilization - general dressing of wounds in out-patient conditions, in OR and hospital wards

ref.	dimensions	grammage	packaging unit/ bag	bulk packaging
NL50-100	5 cm x 5 cm	30 g/m ²	100 pcs	50 bags
NL50-100-40	5 cm x 5 cm	40 g/m ²	100 pcs	50 bags
NL75-100	7.5 cm x 7.5 cm	30 g/m ²	100 pcs	40 bags
NL75-100-40	7.5 cm x 7.5 cm	40 g/m ²	100 pcs	40 bags
NL100-100	10 cm x 10 cm	30 g/m ²	100 pcs	40 bags
NL100-100-40	10 cm x 10 cm	40 g/m ²	100 pcs	40 bags
NL1020-100	10 cm x 20 cm	30 g/m ²	100 pcs	20 bags



NONVI lux S

non-woven swab, sterile

features:

- made of hydrophilic medical non-woven material, 30 or 40 g/m²
- 4 plies
- soft and non-fraying

indication:

- general dressing of wounds in out-patient conditions, in OR and hospital wards



ref.	dimensions	grammage	packaging unit/ bag	intermediate packaging/box	bulk packaging
NL052-S	5 cm x 5 cm	30 g/m ²	2 pcs	25 x 2 pcs	60 boxes
NL053-S	5 cm x 5 cm	30 g/m ²	3 pcs	25 x 3 pcs	45 boxes
NL055-S	5 cm x 5 cm	30 g/m ²	5 pcs	25 x 5 pcs	36 boxes
NL0510-S	5 cm x 5 cm	30 g/m ²	10 pcs	10 x 10 pcs	32 boxes
NL052-S-40	5 cm x 5 cm	40 g/m ²	2 pcs	25 x 2 pcs	60 boxes
NL053-S-40	5 cm x 5 cm	40 g/m ²	3 pcs	25 x 3 pcs	45 boxes
NL055-S-40	5 cm x 5 cm	40 g/m ²	5 pcs	25 x 5 pcs	24 boxes
NL0510-S-40	5 cm x 5 cm	40 g/m ²	10 pcs	10 x 10 pcs	32 boxes
NL752-S	7.5 cm x 7.5 cm	30 g/m ²	2 pcs	25 x 2 pcs	48 boxes
NL753-S	7.5 cm x 7.5 cm	30 g/m ²	3 pcs	25 x 3 pcs	36 boxes
NL755-S	7.5 cm x 7.5 cm	30 g/m ²	5 pcs	25 x 5 pcs	24 boxes
NL7510-S	7.5 cm x 7.5 cm	30 g/m ²	10 pcs	10 x 10 pcs	12 boxes
NL752-S-40	7.5 cm x 7.5 cm	40 g/m ²	2 pcs	25 x 2 pcs	48 boxes
NL753-S-40	7.5 cm x 7.5 cm	40 g/m ²	3 pcs	25 x 3 pcs	24 boxes
NL755-S-40	7.5 cm x 7.5 cm	40 g/m ²	5 pcs	25 x 5 pcs	18 boxes
NL7510-S-40	7.5 cm x 7.5 cm	40 g/m ²	10 pcs	10 x 10 pcs	12 boxes
NL7520-S-40	7.5 cm x 7.5 cm	40 g/m ²	20 pcs	5 x 20 pcs	12 boxes
NL102-S	10 cm x 10 cm	30 g/m ²	2 pcs	25 x 2 pcs	48 boxes
NL103-S	10 cm x 10 cm	30 g/m ²	3 pcs	25 x 3 pcs	24 boxes
NL105-S	10 cm x 10 cm	30 g/m ²	5 pcs	20 x 5 pcs	12 boxes
NL1010-S	10 cm x 10 cm	30 g/m ²	10 pcs	10 x 10 pcs	12 boxes
NL102-S-40	10 cm x 10 cm	40 g/m ²	2 pcs	25 x 2 pcs	30 boxes
NL103-S-40	10 cm x 10 cm	40 g/m ²	3 pcs	25 x 3 pcs	18 boxes
NL105-S-40	10 cm x 10 cm	40 g/m ²	5 pcs	20 x 5 pcs	12 boxes
NL1010-S-40	10 cm x 10 cm	40 g/m ²	10 pcs	10 x 10 pcs	12 boxes
NL1020-S-40	10 cm x 10 cm	40 g/m ²	20 pcs	5 x 20 pcs	12 boxes

Sterilization method:

-S - EO sterilized swabs -SS - steam sterilized swabs